

Client:	Neuro Peptides
Sample:	GLOW Blend 70mg
Batch:	202603GL70

Tests Requested:	Standard + Endotoxins + Sterility
Appearance:	Blue Lyophilized Powder in 3mL Vial

Analysis Conducted on:	03/03/2026	Report Issued on:	08/20/2026
-------------------------------	------------	--------------------------	------------

Test	Method	Result	Specification	
Identity	HPLC-DAD	BPC-157, TB-500, GHK-Cu	BPC-157, TB-500, GHK-Cu	Pass
Purity	HPLC-DAD	99.65	>99% by HPLC	Pass
Quantity	HPLC Analyte Concentration Curve Assay	9.89, 10.54, 57.96mg	10mg, 10mg, 50mg, Respectively	Pass
Endotoxins	LAL ELISA	0.368 EU/mL	≤ 0.5 EU/mL	Pass
Sterility	Pour Plate	No Growth	No Growth	Pass

