

Client: Neuro Peptides**Sample:** Semax 10mg**Batch:** 202603SX10**Tests Requested:** Standard + Endotoxins + Sterility**Appearance:** White Lyophilized Powder in 3mL Vial**Analysis Conducted on:** 03/28/2026**Report Issued on:** 08/20/2026

Test	Method	Result	Specification	
Identity	HPLC-DAD	Semax	Semax	Pass
Purity	HPLC-DAD	99.69	>99% by HPLC	Pass
Quantity	HPLC Analyte Concentration Curve Assay	10.11	10mg	Pass
Endotoxins	LAL ELISA	0.412 EU/mL	≤ 0.5 EU/mL	Pass
Sterility	Pour Plate	No Growth	No Growth	Pass

