

Client: Neuro Peptides

Sample: GHK-Cu 100mg

Batch: 202601GHK100

Tests Requested: Standard + Endotoxins + Sterility

Appearance: Blue Lyophilized Powder in 3mL Vial

Analysis Conducted on: 01/02/2026 **Report Issued on:** 08/20/2026

Test	Method	Result	Specification	
Identity	HPLC-DAD	GHK-Cu	GHK-Cu	Pass
Purity	HPLC-DAD	99.87	>99% by HPLC	Pass
Quantity	HPLC Analyte Concentration Curve Assay	114.54mg	100mg	Pass
Endotoxins	LAL ELISA	0.189 EU/mL	≤ 0.5 EU/mL	Pass
Sterility	Pour Plate	No Growth	No Growth	Pass

