

Client:	Neuro Peptides
Sample:	Methylcobalamin 10mL
Batch:	N/A

Tests Requested:	Standard + Endotoxins + Sterility
Appearance:	Red Liquid in 10mL Vial

Analysis Conducted on:	12/19/2025	Report Issued on:	08/20/2026
-------------------------------	------------	--------------------------	------------

Test	Method	Result	Specification	
Identity	HPLC-DAD	Methylcobalamin	Methylcobalamin	Pass
Purity	HPLC-DAD	N/A	N/A	Pass
Quantity	HPLC Analyte Concentration Curve Assay	.98mg/mL	1mg/mL	Pass
Endotoxins	LAL ELISA	0.124 EU/mL	≤ 0.5 EU/mL	Pass
Sterility	Pour Plate	No Growth	No Growth	Pass

