

**Client:** Neuro Peptides

**Sample:** Semax - 10mg

**Batch:** 202607SX10

**Tests Requested:** Standard + Endotoxins + Sterility Heavy Metals

**Appearance:** White Lyophilized Powder in 3mL Vial


**Analysis Conducted on:** 07/11-12/2026

**Report Issued on:** 08/20/2026

| Test         | Method                                       | Result     | Specification |      |
|--------------|----------------------------------------------|------------|---------------|------|
| Identity     | HPLC-DAD                                     | Semax      | Semax         | Pass |
| Purity       | HPLC-DAD                                     | 99.32%     | >99% by HPLC  | Pass |
| Quantity     | HPLC Analyte<br>Concentration Curve<br>Assay | 12.1mg     | 10mg          | Pass |
| Heavy Metals | Chromogenic Assay                            | ND         | ND            | Pass |
| Endotoxins   | LAL Chromogenic<br>Assay                     | .121 EU/mL | <0.5 EU/mL    | Pass |
| Sterility    | Pour Plate                                   | No Growth  | No Growth     | Pass |



## Chromatogram:



**Disclaimer:** The data contained in this report is representative of the received sample. This analysis is for informational purposes only. The tested substances are intended for research use only and are not approved for human or veterinary use, diagnostic, therapeutic, or clinical applications. All data was obtained under standard laboratory conditions.

